

道  
TAO

**HANGZHOU TAO IMPORT & EXPORT CO.,LTD**



道  
TAO Company sales tenet "Human health, Tao dedication"

道  
TAO HANGZHOU TAO IMPORT & EXPORT CO.,LTD

**Company Profile**



Chairman: Frank Tao

Hang zhou Tao Import & Export Co.,LTD.Registered in beautiful city of Hangzhou,Company production base is located in suburban Tonglu,has about 500 employees. Company products involving medical Endoscope,Lacteal peaceful, Thoracic, Gynecology,Urology,Ent,Orghopedicsetc. Varieties

and specifications,novelstyle,quality is reliable. And for the new old customer customize various special usage and maintenanceof medical instrument and the production of Endos-copic and surgical instruments. To provide us ers with fine quality and reliable service as always, the logistics system.

The company's products in the various provinces, municipalities and autonomous prefectures andclinical trial by experts,won acclaim,believed to be comparable with the same products abroad.

Company sales tenet"Human health,Tao dedication", "Customer need is the direction of our efforts!Off Household Satisfactionis our goal!"Is our belief.



Foreign Trade Director : Manish Kapadia

**Instruments**

<b>Laparoscope Instruments</b>	
Instrument to unload .....	1
Laparoscope .....	1
Clip applicator, TI clip .....	2
Hem-o-lok .....	2
Needle-holding forceps .....	3
Double-electrode forceps .....	4
Electrode instruments. Baffle. Pull. Sping grasping forceps.....	5
Suction&Irrigation tube. Bile duct knife. neele .....	6
Spreader .....	7
Three joint instruments .....	7
Tow joint instruments .....	8
Other instruments .....	9
Scissors,Dissecting forceps .....	10
Other forceps .....	11
Other forceps .....	12
Other forceps .....	13
Trocar,Reducer .....	14
Single-hole(After umbilical) .....	15
Single-hole(After umbilical) .....	15
Single-hole(After umbilical) .....	16
Single-hole(After umbilical) .....	16
Miniature instruments for Laparoscopic surgery .....	17
Instruments for small Incision operation .....	17
Instruments for Thoracoscopic surgery .....	18
Thyroid operating instruments .....	18
Instruments of Biliary tract fetch .....	18
<b>Gynaecology</b>	
Electric uterus-cutter .....	19
Instruments of Gynaecology Laparoscope .....	20
Uterus-cutter (process vagina) .....	21
Hysteroscope .....	21
Cavitas inspectoscope .....	21
Hysterp Resectoscope .....	22
<b>Instrumts for Urological Department</b>	
Cysto-Urethroscopce .....	22
Lithotriposcope .....	23
Urethrotome Telescope .....	23
Instrument of foreskin ring ligator .....	23
<b>Intruments for E.N.T Department</b>	
Sinoscope .....	24
Nose electric-drill .....	24
Hold Laryngoscope .....	25
Esophagoscope instruments .....	25
Otoscope instruments .....	25
Bronchoscope .....	26
<b>Intruments for Orthopaedics</b>	
Arthroscopce instruments .....	26
Endoscope Discectomy System .....	26
<b>Intestine department</b>	
CJ-1 Sigmoidoscope .....	27
CJ-2 Rectoscope .....	27
CJ-3 Sigmoidoscope .....	27

**Other equipments**

Endoscope display systems	28
cold light source	28
insufflator	28
high frequency electric-cutter	28
Irrigation pump hysteroscope	28
E.N., T therapeutic equipments	28
Simulate n. Maneuvering	28
Suction & Irrigation apparatus	29
pUMP	29
Cleaning device for instruments	29
Suction cutter	29

**Operation table**

**Obstetric table**

CC-A Multifunction obstetric table	30
3004 Multi-purpose obstetric table	31
78A Multi-purpose obstetric table	31
DS-8 Electric multi-purpose obstetric operation table	32
DS-3004(check type)Electric multi-purpose obstetric table	33
Model DS-3004 Electric multi-purpose obstetric table	33
Multi-purpose obstetric table	34
Fs common obstetric table	35
<b>Operation table</b>	
DT-2A(translational type)Multi-purpose operation table	36
DT-12C Electric operating table	37
DT-12E Electric operating table	38
DT-12E Electric operating table	39
DT-12F Electric operating table	40
Luxury multifunctional operating table	41
3001 Side operating universal table	42
Head operating universal table	43
JT-SA, JT-2A(translational type)Multi-purpose operation table	44
PT Ordinary operation table	45

**Sterilizer equipment**

Ethylene oxide sterilizer	46
Brief information&Products	47
Products	48
Products	49


**Other**

Wound care products	50
External fixation products	50
Operating room products	51
New products	51


**Instrument to unload**





**Laparoscopes**





Clip applicator.TI Clip




Hem-o-lok




Needle-holding Forceps




Double-electrode forceps



Electrode instruments



Baffle, Pull




Spring grasping forceps



Section & Irrigation tube




Bile duct knife Bile Duct Needle




Spreader




Three joint instruments



Two joint instruments



Other instruments



Scissors Dissecting forceps



101.010  
Dissecting forceps (large)  
Φ 10X330



101.011A  
Dissecting forceps (90°)  
Φ 10X330



101.011B  
Dissecting forceps (90°)  
Φ 5X330



101.012  
Dissecting forceps (straight)  
Φ 5X330



101.013A  
Dissecting forceps (curved)  
Φ 5X330



101.013B  
Dissecting forceps (curved)  
Φ 5X330



101.013C  
Dissecting forceps (curved)  
Φ 5X330



101.014  
Scissors (large curved)  
Φ 10X330



101.016  
Scissors (single-action straight)  
Φ 5X330



101.017  
Scissors (double-action straight)  
Φ 5X330



101.018  
Scissors (single-action curved)  
Φ 5X330



101.019  
Scissors (double-action curved)  
Φ 5X330



101.022  
Hook scissors  
Φ 5X330



101.035  
Grasping forceps  
Φ 5X330



101.036  
Duck-bill forceps  
Φ 5X330



101.036B  
Isthmus grasping forceps  
Φ 5X330




101.037A  
Mid-hollow grasping forceps  
Φ 5X330



101.037B  
Mid-hollow grasping forceps  
Φ 5X330



101.040A  
Bowel forceps  
Φ 5X330



101.040B  
Bowel forceps  
Φ 5X330



101.040B2  
Bowel forceps  
Φ 5X330



101.040C  
Bowel forceps  
Φ 5X330



101.040D  
Bowel forceps  
Φ 5X330



101.041A  
Tumour grasper with teeth  
Φ 5X330

Other forceps



101.041B  
Tumour grasper with teeth  
Φ 5X330



101.041C  
Tumour grasper with teeth  
Φ 5X330



101.042A  
Oviduct forceps  
Φ 5X330



101.042B  
Oviduct forceps  
Φ 5X330



101.042C  
Oviduct forceps  
Φ 5X330



101.043A  
Claw forceps  
Φ 5X330



101.043B  
Claw forceps  
Φ 5X330



101.044  
Grasping forceps (blunt)  
Φ 5X330



101.045  
Cholangiography forceps  
Φ 5X330



101.046  
Parish grasping forceps  
Φ 5X330



101.047  
Retraining tooth forceps  
Φ 5X330



101.047  
Retraining tooth forceps  
Φ 5X330



101.048  
Appendiceal forceps  
Φ 5X330



101.052  
S-knot forceps  
Φ 5X330



101.053  
B-knot forceps  
Φ 5X330



101.054  
Dissecting-knot forceps (curved)  
Φ 5X330



101.064  
Biopsy forceps (with needle)  
Φ 5X330



101.065  
Biopsy forceps  
Φ 5X330




101.076B  
Lung lobe forceps  
Φ 5X330



101.077  
Absorbent cotton forceps  
Φ 5X330



101.078  
Mesentery forceps  
Φ 5X330



101.078B  
Mesentery forceps  
Φ 5X330



101.066B  
Fan-shaped Forceps  
Φ 5X330



101.101  
Fixed grasping forceps  
Φ 5X330

Other forceps



101.101B  
Fixed grasping forceps  
Φ5X330



101.101C  
Fixed grasping forceps  
Φ5X330



101.102A  
Stomach grasping forceps  
Φ5X330



101.102B  
Stomach grasping forceps  
Φ5X330



101.103  
Non-destructive grasping  
forceps  
Φ5X330



101.103B  
Non-destructive grasping  
forceps  
Φ5X330



101.104  
Octopus grasping forceps  
Φ5X330



101.105  
hook forceps  
Φ5X330



Other forceps

Φ10mm



101.005B  
Biopsy forceps Φ10x330



101.006  
Grasping forceps Φ10X330



101.006B  
Grasping forceps Φ10x330



104.003  
Fastness forceps(thick tooth)  
Φ10X330



104.004  
Broadness mouth forceps  
Φ10X330



104.006  
Atraumatic forceps  
Φ10X330



104.007  
Grasping forceps (long-jaw)  
Φ10X330



104.009  
Duck-jaw forceps  
Φ10X330



101.007A  
Intestinal forceps Φ10x330



101.076A  
Lung lobe forceps Φ10X330



101.079  
Target forceps  
Φ10X330




101.090  
Suture forceps Φ10x330



101.007B  
Intestinal forceps Φ10x330


Trocar, Reducer




KE-101DK Single-hole (After umbilical)



KE-101DK Single-hole (After umbilical)



KE-101DK Single-hole(After umbilical)



KE-201C



KE-102 Miniature instruments for Laparoscopic surgery



KE-103 Instruments for small Incision operation



KE-104 Instruments for thoracoscopic surgery

CODE	NAME	REFURON	QUANTITY
101.001	Clip applicator (mini)	Ø10X330	1
101.076A	Lung lobe forceps	Ø10X330	1
104.003	Full-thickness forceps(black tooth)	Ø10X330	1
104.004	Broad-based mouth forceps	Ø10X330	1
101.077A	Panhandle fan-shaped forceps	Ø10X330	1
104.006	Atraumatic forceps	Ø10X330	1
104.007	Grasping forceps (long-jaw)	Ø10X330	1
101.006	Biopsy forceps	Ø10X330	1
104.008	Double-jaw forceps	Ø10X330	1
201.021	Blot duct needle	Ø5X330	1
102.020B	Knit Pusher	Ø5X330	1
101.043A	Claw forceps	Ø5X330	1
101.043B	Claw forceps	Ø5X330	1
101.017	Scissors(double-action straight)	Ø5X330	1
101.018	Scissors(double-action curved)	Ø5X330	1
101.022	Hook scissors	Ø5X330	1
101.012	Dissecting forceps(straight)	Ø5X330	1
101.013	Dissecting forceps(curved)	Ø5X330	1
101.018B	Dissecting forceps(90°)	Ø5X330	1
101.027	Hook- electrode	Ø5X330	1
101.028	Bipolar-electrode	Ø5X330	1
101.029	Knife- electrode	Ø5X330	1
101.028	Spatial-electrode	Ø5X330	1
101.031	Needle-electrode	Ø5X330	1
104.028	Trocar	Ø8.5	1
104.028	Trocar	Ø8.5	1
104.028	Trocar	Ø8.5	1
104.027	Trocar	Ø12.5	1
104.031A	Irigation tube	Ø10X330	1
104.031B	Irigation tube	Ø5X330	1
104.032	Lymph forceps	Ø5X330	1
101.027	Lymph forceps	Ø5X330	1

KE-106 Thyroid operating instruments

CODE	NAME	REFURON	QUANTITY
101.019	Scissors(double-action curved)	Ø5X330	1
101.013	Dissecting forceps(curved)	Ø5X330	1
106.003	Trocar	Ø10.5X125	1
106.004	Trocar	Ø5.5X125	1
101.071A	Retractor	Ø10.0-5	1
101.026	Needle-holder forceps(retractor)	Ø5X330	1
106.007	Stripper	Ø5X400	1
106.008	Injection needle	Ø2.5X400	1
106.009	Thyroid dissecting	Ø10X400	1
106.010	Stripper(curved)	Ø5X400	1
106.011	Hook	Ø5X330	2
101.068	suction & irrigation tube	Ø5X330	1

KE-107 Instruments of Biliary tract fetch

CODE	NAME	REFURON	QUANTITY
107.001	Vain forceps	Ø10X330	1
107.002	Vain inspection forceps	Ø10X330	1
107.003	Biliary tract fetch forceps	Ø5X330	1
107.004	Trocar(needle)	Ø12X95	1

KE-201A Electric uterus-cutter

CODE	NAME	REFURON	QUANTITY
201.001	Controller		1
201.002	Motor with Handpiece		1
201.003	Grasping forceps(large)	Ø10X420	1
201.004	Grasping forceps(medium)	Ø8X420	1
201.005	Trocar sleeve	Ø18X96	1
201.006	Trocar sleeve	Ø10-Ø18	1
201.007	Trocar sleeve	Ø10-Ø15	1
201.008	Dilator	Ø10X330	1
201.009	Guide probe	Ø10-Ø10	1
201.010	Retractor	Ø18/Ø10	1
201.011	Retractor	Ø18/Ø10	1
201.012	Morcellator blades	Ø18X330	1
201.013	Morcellator blades	Ø15X330	1
201.014	Morcellator blades(cavity)	Ø18X330	1
201.015	Morcellator blades(cavity)	Ø15X330	1
201.016	Ligator	Ø5X330	1
201.017	Hysteromyoma bodkin	Ø10X450	1
201.018	Womb measuring probe	Ø5X500	1
201.019	Palpation probe	Ø10X400	1
201.020	Cervix forceps(adjustable)	250	1
201.021	Transfer needle	Ø5X330	1
KE-002	Sealing cap		10
201.022	Morcellator blades	Ø10X330	1
201.023	Morcellator blades(cavity)	Ø10X330	1
201.024	Apparatus	400	1
201.025	Cervix forceps	Ø5X420	1
201.027	Hysteromyoma bodkin	Ø5X420	1
201.027	Palpation probe	Ø5X420	1

KE-201C Instrument of Gynaecology Laparoscope

201.043 Mightiness-Apparatus




201.044 Feminine-Apparatus



201.045 Curved-Apparatus



201.046 Yaogens-Apparatus



KE-201B Uterus-cutter(process vagina)

201.028  
Vaginal depressor



201.030  
Myoma stripping




201.031  
Ligament hook-shape forceps




201.032  
Single grasping forceps



201.033  
Double grasping forceps



KE-202 Hysteroscope



CODE	NAME	OPERATOR	QUANTITY
202.001	Endoscope (s27)	Φ4X310	1
202.002	Outer sheath	FD24	1
202.003	Working element		1
202.004	Chiselator		1
202.005	Blopsy forceps		1
202.006	Scissors		1
202.007	Grasping forceps		1
202.008	Foreign body forceps		1
202.009	Brush		1
301.014	Intake connector		1
301.015	Outside connector		1
KE-001A	Fiber optic light cable	4X1800	1
KE-604	Pump		1

KE-203 Cavitas inspectoscope

203.001 22' Cavitas inspectoscope



KE-204 Hystero Resectoscope

CODE	NAME	PERFORM	QUANTITY
204.001	Endoscope	Φ4X310	1
204.002	Inner sheath	3600旋轉	1
204.003	Outer sheath		1
204.004	Obturator		1
204.005	Working element	鑷形式	1
204.006	Cutting loop		4
204.007	Roller electrode		2
204.008	Ring electrode		4
204.009	Distraction-roller electrode	500W	2
204.010	Hyperator		1
301.014	Intake connector		1
301.015	Outside connector		1
KE-001A	Fiber optic light cable		1

KE-301 Cysto-Urethroscope

CODE	NAME	PERFORM	QUANTITY
301.001	Endoscope 0°	Φ4X302	
301.002	Endoscope 30°	Φ4X302	
301.003	Endoscope 70°	Φ4X302	
301.004	Working element		1
301.005	Sheath	Φ18.5X230	
301.006	Sheath	Φ19.2X230	
301.007	Sheath	Φ25X230	
301.008	Endoscope bridge		
301.010	Flexible biopsy forceps	Φ2.3X400	
301.011	Flexible grasping forceps	Φ2.3X400	
301.012	Flexible biopsy forceps	Φ2.3X400	
301.013	Flexible grasping forceps	Φ2.3X400	
301.014	Intake connector		
301.015	Outside connector		
301.016	Cleaning stick		
301.017	Fiber optic light cable	Φ4X1800	
KE-002	Sealing cap		

KE-304 Lithotripsy

CODE	NAME	PERFORM
304.001	Endoscope 70°	Φ4
304.002	Lithoclast	
301.014	Intake connector	
301.015	Outside connector	Φ18.5X230
304.004	Sheet	
304.005	Intigating connector	
KE-001	Fiber optic light cable	Φ4X1800


KE-305 Urethrotome Telescope

CODE	NAME	PERFORM
305.001	Endoscope 12°	Φ4x302
305.002	Working handle	Passive
305.003	Sheath	
305.004	Adapter	
305.005	Lance blade	
305.006	Half-moon blade	
305.007	Sickle-shaped blade	
301.014	Intake connector	
301.015	Outside connector	
KE-001A	Fiber optic light cable	Φ4x1800

KE-037 Instrument of foreskin ring ligator

CODE	NAME
307.001	Diagonal rolling forceps
307.002	Inside diameter forceps
307.003	Outside diameter forceps
307.004	scissors
307.005	fat cutting forceps
307.006	conduct out-rlng forceps
307.007	Ring


KE-501 Sinoscope




KE-503 Nose electric-drill



KE-502 Hold Laryngoscope




KE-506 Esophagoscope instruments




KE-505 Otoscope Instruments




KE-504 Pediatric Bronchoscopy sets



KE-401 Arthroscope



KE-402 Endoscopic Discectomy system



Sigmoidoscope



KE-601 Endoscope display systems



Other equipment



KE-603 Suction & Irrigation apparatus



Suction cutter

X100-1

Suction cutter JRX-I

X100-2

Suction cutter JRX-I



## CC-A MULTIFUNCTION OBSTETRIC TABLE

### CC-A Type CC-A

Both the body table and leg table are movable lockable.



CC-A multifunctional obstetric table is designed for gynecological inspection, diagnosis and operation, and is a special chair for delivery as well. The board width makes possible to accommodate both mother and newborn baby so that a good nursing can be guaranteed. The movement of this chair is controlled manually. The new imported gas spring controls positions of body table, so it is easy to operate, convenient and safe. Lifting of body table and foot wheels is conducted through a hydraulic system and transversal valve, which could lessen work intensity of medical staff. Adjustment of leg table is driven by mechanical gearing. It is equipped with a lot of accessories, so it is a multi-functional obstetric table.

#### Main specification parameters

Item	Contents
Length of body table and foot table	1150mmX850mm
Adjustment range of back section	up 0-42° down 0-6°
Adjustment range of abdomen section	up 0-20° down 0-8°
Max. degree of forward and backward tilting of table	forward ≥8° backward ≥22°
Lowest height of body bed	640mm-880mm
Size of table board	body table 1150X650mm leg table 850X600mm

## 3004 MULTI-PURPOSE OBSTETRIC TABLE



#### Main specification parameters

Item	Contents
Length and width of table	1800mm x 600mm
Height adjustment range	Min 650mm x Max 895mm
Sizes of sections	
Back section	600mm x 800mm
Seat section	600mm x 400mm
Leg section	540mm x 500mm
Max. degree of back rest adjustment	up >75° down >6°
Max. degree of forward and backward tilting of tabletop	forward >10° backward >25°

## 78A MULTI-PURPOSE OBSTETRIC TABLE



#### Main specification parameters

Item	Contents
Length and width of table	2000mmX600mm
Height adjustment range	Min 700mm Max 940mm
Sizes of sections	
Back section	850mmX600mm
Seat section	400mmX600mm
Leg section	750mmX600mm
Max. degree of back rest adjustment	up >80° down >8°
Max lifting height of waist board	120mm
Height adjustment range	Min 700mm Max 940mm

78A obstetric table is suitable for delivery and various gynaecological surgeries. It has advantages such as multi-function, convenient manipulation, safe and reliability. The whole table is composed of the main bed and table-leg bed. Lifting of the main bed is controlled by hydraulic pump, and table-leg bed is gas spring.

DS-8 ELECTRIC MULTI-PURPOSE OBSTETRIC OPERATION TABLE

DS-8 / Type DS-8



DS-8 electric multi-purpose obstetric operation table is a kind of comprehensive operating table for delivery, gynecological surgery, diagnosis and examination. The table is equipped with mini-touch controllers to operate lifting, front and back inclining, upward and downward folding of back rest, and fixing and moving of base. It is easy to operate and convenient. Its key parts used imported electric machine with low noise and stable performance.  
Adjustment of leg table is driven by mechanical gearing. It is equipped with a lot of accessories, so it is a multi-functional obstetric table.

Main specification parameters

Item	Contents
Length and width of table	Main bed 1150mmX650mm leg table 850mmX600mm
Max. & min. Height of tabletop	900mm x600mm
Max. degree of forward and backward tilting of table	forward $\geq 10^\circ$ backward $\geq 25^\circ$
Max. degree of back rest adjustment	up $\geq 75^\circ$ down $\geq 15^\circ$
Power supply	~220V 50Hz 1.0Kw

DS-3004(CHECK TYPE)ELECTRIC MULTI-PURPOSE OBSTETRIC TABLE

DS-3004 / DS-3004(Check type)



DS-3004 (Check type) electric multi-purpose obstetric table is powered by electric lead screws. It is very convenient to control the up-and-down movement of table and the adjustment of the tilt angle between back rest and seat by foot/manual switch.

If want to make check the bed move, it is only to loosen slightly the tight device of bottom four feet and tight slightly after moving to the needed position. Foot mattress is an additional device which can be put on or off conveniently according to the surgical needs.

Main specification parameters

Item	Contents
Length and width of table	1900mm x 600mm
Height adjustment range	Max $\geq 900$ mm Min $\leq 500$ mm
Max. degree of back rest adjustment	-15° ~70°
The maximal rotation angle of dicky	-10° ~25°
Power supply	220V 50HZ 110VA

MODEL DS-3004 ELECTRIC MULTI-PURPOSE OBSTETRIC TABLE

DS-3004 / multi-purpose obstetric table



DS-3004 multi-purpose obstetric table is designed for delivery, abortion and gynecological examination. Lifting of the bed and back rest and front and back inclining of table are controlled by advanced electric push rod and foot switch. The fix-movements of base is controlled by foot pedal.

Main specification parameters

Item	Contents
Length and width of table	1900mm x 600mm
Height adjustment range	Min 650mm x Max 895mm
Sizes of sections	
Back section	600mm x 800mm
Seat section	600mm x 400mm
Leg section	600mmX600mm
Max. degree of back rest adjustment	up $>75^\circ$ down $>6^\circ$

MULTI-PURPOSE OBSTETRIC TABLE

DS-3004A



3004



FS COMMON OBSTETRIC TABLE



FS

The bed is used for delivery and various gynecological surgeries, as well as family planning. It can be used for parturients' delivery and surgery with various postures by appropriate operation. One of its characteristics is multi-function. It has adopted separation form of main bed and leg bed in the structure. So when moving the parts of leg bed, it is greatly convenient for doctors to operate surgery. Back section and seat section are controlled by hand wheel and handle. It is easy to operate, reliable and stable, so it is the ideal, necessary gynecological instrument for various hospitals.

Main specification parameters

Item	Contents
Length and width of table	1950mmX600X800mm( Height )
Sizes of sections	
Back section	750mmX600mm
Seat section	400mmX600mm
Leg section	750mmX600mm
Max. degree of back rest adjustment	up $\geq 60^\circ$ down $\geq 8^\circ$
Max. degree of seat section adjustment	up $\geq 25^\circ$ down $\geq 8^\circ$



FS-B Stainless Steel Gynecologic Examination Bed

Main specification parameters

Item	Contents
Length and width of table	1400mmX800mmX570mm( Width )
Max. degree of back rest adjustment	down $\geq 18^\circ$



DT-2A(TRANSLATIONAL TYPE)MULTI-PURPOSE OPERATION TABLE



DT-2A

Optional forked leg plate  
Power voltage, power frequency,  
power capacity  
~ 220V , 50Hz , 1.0KW

- The controller device lock key protecting from misoperation.
- The Table is made of penetrated X-ray material. There is cassette guide rail under the table which make radiography.
- The table is fixed a button which can control the electric wheel to draw in and put out so that the table can be pulled easily.
- Accessories and Guide rail are made of stainless steel, polyurethane etc material.
- The main accessory hydraulic pump and electromagnetic valve are made of imported component from Taiwan.
- The mattress is made of foam polyurethane , acid and alkali resistance , high temperature resistance and easy clean.
- The head board and leg board are adopted special imported gas-spring from Germany . it is convenient to operate.

Main specification parameters

Length and width of table	2000mm x 480mm
Max. &Min. Height of table	930mm x 690mm
Max. degree of forward and backward tilting of table	forward $\geq 25^\circ$ backward $\geq 20^\circ$
Max. degree of left and right tilting of table	left $\geq 18^\circ$ right $\geq 18^\circ$
Range of leg rest adjustment	down $\geq 90^\circ$ detachable
Range of back rest adjustment	up $\geq 75^\circ$ down $\geq 10^\circ$
Range of head rest adjustment	up $\geq 45^\circ$ down $\geq 50^\circ$ detachable
Elevating height of kidney bridge	$\geq 80$ mm

DT-2A



DT-12C ELECTRIC OPERATING TABLE



- The bedstead is made of aluminium alloy.
- The accessories , base shell and column shell all are made of stainless steel.
- The main accessory hydraulic pump and electromagnetic valve are adopted imported component
- The tabletop is made of Paper plywood (double-layer), there is a film-frame in the interlayer, it is convenient to Move the radiography Position.
- The leg board can be forked and abducted  $180^\circ$  , folded down  $90^\circ$  . It is suitable for urinary 、 obstetrics and Gynecology surgery

Main specification parameters

Item	Contents
Length and width of table	2050mmX480mm
Max. &Min. Height of table	950mmX700mm
Max. degree of forward and backward tilting of table	forward $\geq 20^\circ$ backward $\geq 20^\circ$
Max. degree of left and right tilting of table	left $\geq 18^\circ$ right $\geq 18^\circ$
Range of back rest adjustment	up $\geq 75^\circ$ down $\geq 20^\circ$
Max. degree of leg rest adjustment	$\geq 90^\circ$ abduction $180^\circ$
Range of head rest adjustment	up $\geq 75^\circ$ down $\geq 90^\circ$ detachable
waist board	$\geq 120$ mm
power	~220V 50Hz 1.0Kw

DT-12E ELECTRIC OPERATING TABLE



Right manipulation

- The table is equipped with microcomputer controller and safety lock key protecting from misoperation.
- Adjusted manipulably to the enough low place (the lowest height 550mm) so that doctor can sit to operate cerebral surgery ophthalmology.
- The Table can be moved longitudinally more than 300mm so that C arm X-ray can rotate freely to work ,and x-ray can make radiography vertically without dead angles.
- The Table is divided into five section : head board , shoulder board、back board 、 seat board and leg board.
- The Table is made of penetrated X-ray material. There is cassette guide rail under the table which make radiography.
- Equipped with shoulder-back composited bending button on the table , which will be convenient for kidney surgery and Gallbladder surgery.
- The table is fixed a button which can control the electric wheel to draw in and put out so that the table can be pulled easily.
- Accessories and Guide rail are made of stainless steel.
- The main accessory hydraulic pump and electromagnetic valve are made of imported component from Taiwan.
- Base shell is made of aluminum alloy and supporting base is made of stainless steel.
- The head board and leg board is adopted special gas-spring from Germany, it is convenient to operate.

DT-12E ELECTRIC OPERATING TABLE



Touch-sensitive display

Type DT-12E

It can be equipped with two independent electric hydraulic systems

- It adopts imported electric appliances and electric hydraulic system.
- The table is constructed with 5 movable parts, and suitable for all kinds of surgeries as well as special operative positions, and C arm X-ray inspection and radiography.
- It is equipped with X-ray rack and cassette specifically for C arm users.
- It is equipped with microcomputer controller and safety lock key protecting from mishandling. Set one-key restoration device.
- Unilateral column design and large space under table have given adequate space to install C arm X-ray apparatus and other device.
- Steering control switch and universal steering wheel make the table move and rotate free when needed.
- Removable and assemble head board and leg plate which can be forked and lengthened. Generous, nice looking and humanization design make clinicians stable gravity centre and decrease fatigue stress during operation.
- Shell is made of stainless steel. Mattress is made using new technology and new materials, so it has many advantages such as easy clean, acid and alkali resistance, anti-bacterial, impact, comfortable and beautiful.

DT-12F ELECTRIC OPERATING TABLE



- The mattress is made of foam polyurethane, acid and alkali resistance, high temperature resistance and easy clean.
- The double-layer tabletop is suit for X-ray examination when you operate. The leg board can be abducted 180°, folded 90° and detached.
- The wheels adjustment, table height adjustment, forward and backward tilting left and right tilting, back board folded are all controlled by linear motor. It's convenient and safe.
- Accessories, guide rail, base shell and column shell all are made of stainless steel. The table is equipped Auto-Charge device; It can keep the table being used when emergency power shuts off.

DT-12F



LUXURY MULTIFUNCTIONAL OPERATING TABLE

1088



3008H

3008H LUXURY MULTIFUNCTIONAL OPERATING TABLE

- The bedstead, accessories, guide rail, base shell and column shell all are made of stainless steel, and antistatic finish as well.
- Back board and leg board are controlled by imported gas spring from Germany.
- The mattress is made of foam polyurethane, it has many advantages such as Easy clean, acid and alkali resistance, high-Temperature resistance



### 3001 SIDE OPERATING UNIVERSAL TABLE

The operating table is supplied for thoracic surgery, abdominal surgery, ophthalmology, otorhinolaryngology, obstetrics and gynecology, urinary surgery and orthopaedics.

Lifting of this operating table is controlled by foot oil pump, and the rest of the movements are controlled by mechanical transmission side manipulation.

#### Main specification parameters

Length and width of table	2100mm×480mm
Max.&Min.Height of table	1045mm×800mm
Max.degree of forward and backward tilting of table	forward ≥35° backward ≥20°
Max.degree of left and right tilting of table	left ≥22° right ≥22°
Oil pump travel	240mm
Head board (275×310mm)	It can be folded for 90 degree, elongated or unloaded.
Back board (500×480mm)	down ≥22° , up 75°
Kidney bridge (85×415mm)	It can be increased up to 120mm, or equal to the table surface.
Leg plate (500×480mm)	down ≥90°
Seat section (450×480mm)	There is a round hole between seat rest and leg rest.



● TYPE3001 ●



#### 3001-1



### HEAD OPERATING UNIVERSAL TABLE

The operating table is supplied for thoracic surgery, abdominal surgery, ophthalmology, otorhinolaryngology, obstetrics and gynecology, urinary surgery and orthopaedics.

The rise and drop of operating table was handled by pedaling the oil pump, other various movements were handled by the mechanical transmission head.



●3008D fored leg plate



●3008TYPE●



3008A TYPE

Base shell made of stainless steel  
Additional film-frame

3008B TYPE

Base, shell and bed surface made of stainless steel  
Additional film-frame

#### Main specification parameters

Item	Contents
Length and width of table	2100mm×480mm
Max.&Min.Height of table	1045mm×800mm
Max.degree of forward and backward tilting of table	forward ≥55° back ward ≥35°
Max.degree of left and right tilting of table	left ≥25° right ≥25°
Max.degree of back rest adjustment	up ≥85° down ≥10°
Waist rest elevation	120mm
Leg rest adjustment	down ≥90°
Oil pump travel	240mm
Head rest (275mm×310mm)	It can be folded for 90 degree, elongated or unloaded.

JT-2A JT-2A(TRANSLATIONAL TYPE)MULTI-PURPOSE OPERATION TABLE



Type JT-2A



Optional forked leg plate

This multi-functional operating table is used for operations and examination of thoracic surgery, abdominal surgery, orthopaedics, cerebral surgery, ophthalmology, otorhinolaryngology, obstetrics and gynecology and urinary surgery. There is a shelf fixed under the table. Cassette is put on the shelf and adjusted to the radiography position, then X-ray radiography is taken. The mattress is made of foam polyurethane. The bottom seat is movable and can be locked. Using the imported air-spring, the table has got the best position for clinical uses with its unique electric hydropneumatic technology and is convenient to operate, so it is helpful to improve operation efficiency of doctors.

Main specification parameters

Item	Contents
Length and width of table	1960mmX480mm
Max.&Min.Height of table	930mmX690mm
Max.degree of forward and backward tilting of table	forward $\geq 25^\circ$ backward $\geq 20^\circ$
Max.degree of left and right tilting of table	left $\geq 18^\circ$ right $\geq 18^\circ$
Range of leg rest adjustment	down $\geq 90^\circ$ detachable
Range of back rest adjustment	up $\geq 75^\circ$ down $\geq 10^\circ$
Range of head rest adjustment	up $\geq 45^\circ$ down $\geq 50^\circ$ detachable
Elevating height of kidney bridge	$\geq 80$ mm

JT-2A



PT ORDINARY OPERATION TABLE



PT ORDINARY OPERATION TABLE



PT Ordinary Operation Table

The operating table is supplied for head surgery, neck surgery, thoracic surgery, abdominal surgery, vulvar surgery and limbs surgery as well as obstetrics and gynecology, ophthalmology, otorhinolaryngology and orthopaedics.

Main specification parameters


Item	Contents
Length and width of table	2000mmX480mm
Height adjustment range	Min $\leq 800$ mm Lifting height <sup>off table</sup> $\geq 250$ mm
Leg rest	down $\geq 90^\circ$
Back rest	up $\geq 90^\circ$ down $\geq 20^\circ$
Kidney bridge	Lifting height $\geq 120$ mm
Head rest	up $\geq 45^\circ$ down $\geq 90^\circ$

The operating table is supplied for head surgery, neck surgery, thoracic surgery, abdominal surgery, vulvar surgery and limbs surgery as well as obstetrics and gynecology, ophthalmology, otorhinolaryngology and orthopaedics.



## EHTYLENE OXIDE STERILIZER

Professional ETO Sterilizer Manufacture



TAO GROUP

TAO Voice Technology

TAO Import & Export Co Ltd (HSDP Group)

HEAD office: Room1602-1604, Fortune Plaza, No.63, North Huan Cheng Road, Hangzhou, 310004, China

TEL: Direct +86-571-28281830 /28281805 Fax: 86-571-28281800

Mobile: 86-18368038308 / 18958035786

email:manish@taoandcompany.com / manish@frankdesigner.com

## BRIEF INFORMATION

We are the leading manufacturer of OXIDE STERILIZER in china , Our Tao Group series sterilizer which is conformity with CE Standard have been widely used in Medical Device , health care Facilities plastics and Disposables products, culture goods and currency etc. We also could manufacture matched Pre-conditioning system, aeration room and professionally ethylene oxide residual gas treatment system.

Our Technical team (Tao Group) is working in Sterilization fields more than 20 years and we are the drafter of China {ETHYLENE OXIDE STERILIZATION} Standard.

Our Sales team {TAO GROUP} has a lot of experience in Pre-Sales and after Sales Service. We have more than 70% Market Shares in china and our products have been exported to many countries like EGYPT .TURKEY,SOUYTH AFRICA , ALGERIA ,PARAHAUY , UZBEKISTAN ,MALAYSIA , RUSSIA ,POLAN ,VENEZUELA,etc .

We are looking forward to establishing friendly business relation under mutual benefit between our companies from now on. Give you our best goods and perfect service is our promises WELCOME TO OUR FACTORY AND SERVICE QUALITY IS OUR MOTTO

## PRODUCTS

We also can manufacture special Size Sterilizer according to your real demands and requirements beyond the following specification

SPECIFICATION / SIZE	OUTSIDESIZE (LXBXH)mm	TRANSPORT SIZE {LxBxH}mm	INSIDE SIZE {LxBxH}mm
TAO-EN-25/25MTR SQ	10450X3050X2300	10450X2000X2300	10450X1450X1700
TAO-EN-20/20MTR SQ	9450X2950X2300	9450X1900X2300	9450X1350X1700
TAO-EN-15/15 MTR SQ	7200X2950X2300	7200X1900X2300	7200X1350X1700
TAO-EN-10/10MTR SQ	4950X2950X2300	4950X1900X2300	4500X1350X1700
TAO-EN-6/6MTR SQ	3250X2950X2300	2300X2300X2300	1850X1000X1650
TAO-EN-1/1MTR SQ	1550X1800X1800	1550X1800X1800	1260X800X1050
TAO-EN-0.2/0.2MTR SQ	1000X850X1500	1000X850X1500	650X550X550

## PRODUCTS

### 2) Working Principle

Ethylene Oxide is a Broad-Spectrum bactericidal agent and can be used to kill various kinds of micro-organism including bacterial Vegetativeform, spore, virus, and fungal spore. It is generally agreed that Ethylene Oxide kills the micro-organism because it produces a non-special alkylation with the protein, DNA, RNA.

### 3) Advantage

EQ Sterilization is a Chemical process consisting of four primary variables, GAS concentration, Humidity, Temperature and time. In comparison to other sterilization.

EQ sterilization is reliability and gentleness ,high permeable, board-spectrum bactericide which can kill all known viruses , bacteria ,and fungi.annihilation even the most Sterilization-resistant types of micro-organism, bacterial spores.

### 4) Main Structure

Ethylene Oxide System is consisting of sterilizing chamber, vacuum system, heating circulation system. Humidification system, gasifying system, door and seal system. Residual gas treatment system etc.

### 5) Primary Parameter of Technical

Sterilizing temperature: 20 – 100 degree C (empty chamber+<sub>3</sub> degree C, Full chamber<10Degree C)

Sterilizing Pressure:-85K Pa - + 85K Pa Precision: 1 Grade, resolution: 0.1KPa

Sterilizing Humidity: 30% RH-85% RH

Humidification: Arbitrarily

Sterilization concentration: 250mg/L, 1200 mg/L

Sterilizing time: 30 min -12 hours

Max power :< 65KM, Electrical Heating

Max Lead rate :< 0.1KPa/MIN

Power Source: 3 phase, 4 wire, 380 V <sub>+10%</sub>, 50Hz

Working temperature: 5-40Degree C, Humidity :< 80%RH



## PRODUCTS

### 6) Applied Range

#### 1. Medical Device:

Injector , infusion equipment, surgical gloves, various catheter, oxygen breathing apparatus ,suture ,operation kits, anesthetic facilities ,medical dressing , endoscope, cardiac pacemaker, dental device, heart-lung machine etc.

#### 2. Health Care Article

Women sanitary Towel, Napkin, Disposable Chopsticks, Disposable Dishware

#### 3. Medical Material and cosmetic

Chinese medicinal plant, capsule, western medicines and partial cosmetic articles etc.

#### 4. Culture Goods and Currency

Bibliographicfiles, books, and paintings, cultural relic, currency and herbariumetc. Clothing, leather, pelage, and fur, spices, chemical fiber, Teryleneetc. . . .



Wound care products



First aid kit



Medical Dressing Kits



Lap Sponges

External fixation products



Tubular Elastic Bandages



Universal Elastic Bandages



Crepe Bandages



Knitted Elastic Bandages



Gauze bandage

Operating room products



Visitor Coat



SMS Coverall



Microporous Coverall Plus Tape Back

New products



varicose veins care arm

Compression Stocking  
Compression



Antithrombotic Stockings



## HANGZHOU TAO IMPORT & EXPORT CO.,LTD

Office Address: Room 1602-1604, Fortune Plaza, No.63 North Huancheng Road, Hangzhou, China

Showroom Address: 10F, Kanghua Building, No.100 GongKang Road, HangZhou, China Postcode: 310012

Tel: 0086 571 85273789 / 0086 571 28281830/805

Mobile: +86 18368038308 / 18958035786

E-mail: [manish@taoandcompany.com](mailto:manish@taoandcompany.com) / [manish@frankdesigner.com](mailto:manish@frankdesigner.com) / [gurongxuan@taoandcompany.com](mailto:gurongxuan@taoandcompany.com)

